



CABLE GLAND SELECTOR TABLE

The gland sizes are only a recommendation.
Actual sizes may vary slightly between different brands of gland.

CABLE SPECIFICATION		GLAND SIZE		
Conductor Area mm ²	Number of Cores	BW Indoor	CW Outdoor	E1W & E1WF
1.5mm	2	20s	20ss	20ss
	3	20s	20ss	20ss
	4	20s	20s	20s
	7	20s	20s	20s
	12	20	20	20
	19	25	25	25
	27	25	25	25
	37	32	32	32
	48	32	32	32
2.5mm	2	20s	20s	20s
	3	20s	20s	20s
	4	20s	20s	20s
	7	20	20	20
	12	25	25	25
	19	25	25	25
	27	32	32	32
	37	40	40	40
4mm	2	20s	20s	20s
	3	20	20	20
	4	20	20	20
6mm	2	20	20	20
	3	20	20	20
	4	20	20	20
10mm	2	25	25	25
	3	25	25	25
	4	25	25	25
16mm	2	25	25	25
	3	25	25	25
	4	25	25	25
25mm	2	25	25	25
	3	32	32	32
	4	32	32	32
50mm	2	32	32	32
	3	32	32	32
	4	40	40	40
70mm	2	32	32	32
	3	40	40	40
	4	40	40	40

95mm	2	40	40	40
	3	40	40	40
	4	50	50	50s
120mm	2	40	40	40
	3	50	50	50s
	4	50	50	50
150mm	2	50	50	50s
	3	50	50	50s
	4	50	50	50
185mm	2	50	50	50s
	3	50	50	50
	4	63	63	63s
240mm	2	50	50	50
	3	63	63	63s
	4	63	63	63
300mm	2	63	63	63s
	3	63	63	63
	4	75	75	75s
400mm	2	63	63	63
	3	75	75	75s
	4	75	75	75

